

# CRANIOCORPOGRAPHY



## Fast and simple investigation


NeuroEquilibrium™ Craniocorpography (CCG) is a neuro-otological investigation to evaluate the vestibulospinal reflex, which is responsible for maintenance of balance during gait testing.

The protocols evaluated by this screening include

- Romberg test
- Tandem walking
- Unterberger's test, which evaluates the following parameters
  - Displacement
  - Sway
  - Angle of deviation
  - Body axis spin

## Product specifications

- Real time data analysis
- Fully integrated and computerized
- Patent pending



## Cloud integration for anywhere access

• Data from the tests is automatically and seamlessly saved on the cloud and can be accessed over internet from anywhere in the world.

• Allows for easy restoration of documentation in case of local system crashes.



NeuroEquilibrium™ provides a comprehensive management of vertigo and balance disorders by breakthrough innovations. Our strategy is to set up specialized advanced vertigo and balance clinics across India, Asia and Africa to treat the huge numbers of vertigo and balance disorder patients through our

- Neurotological evaluation with our suite of world class diagnostic equipment;
- Specialist neurotological team; and
- NeuroEquilibrium™ Assessment Platform making it India's first robust diagnostic protocol to help doctors diagnose and treat vertigo cases easily across continents.

Let's **get together** to  
restore **balance**™



  
**NeuroEquilibrium**

Diagnostic Systems Pvt Ltd

140 Sri Gopal Nagar  
Gopalpura Bypass Road  
Jaipur 302015  
India  
+91- 8003897909  
info@neuroequilibrium.in  
www.neuroequilibrium.in



M4570910IN  
ISO 13485:2003



---

Video Nystagmography  
(VNG)

---

Subjective Visual Vertical  
(SVV)

---

Dynamic Visual Acuity  
(DVA)

---

NeuroEquilibrium™  
Assessment Platform