



Neuro Equilibrium[®]

**Cutting-edge technology for
diagnosis and treatment of
vertigo & balance disorders**

CALORIC IRRIGATOR

Caloric test has been the gold standard in vestibular testing for over 100 years. NeuroEquilibrium Caloric Irrigators are used in combination with the VNG machine to conduct caloric tests on patients with dizziness or balance disorders

Water Caloric Irrigator

provides water at specified flow and fixed temperatures for caloric testing of the patient.

Features

- Ergonomically designed delivery head
- Visual cues indicating beginning and end of irrigation
- Long 9-ft delivery hose to easily reach the patient's ear across an exam table
- External water tank that holds enough water to test three patients
- Water flow can be started with delivery head switch

Product Specifications

- Cold stimulus temperature at 30°C
- Warm stimulus temperature at 44°C
- Temperature accuracy +/- 1°C
- Flow rate: 500 ml per minute
- Computer Integration: Integrated with NeuroEquilibrium VNG



Features

- Ergonomically designed delivery head
- Visual cues indicating beginning and end of irrigation
- Long 9-ft delivery hose to easily reach the patient's ear across an exam table
- Air flow can be started with delivery head switch

Air Caloric Irrigator

provides air at specified flow and fixed temperatures for caloric testing of the patients with vertigo and balance disorders. The Air Caloric Irrigator is especially useful in cases of ear perforation.

Product Specifications

- Cold stimulus temperature at 24°C
- Warm stimulus temperature at 50°C
- Temperature accuracy +/- 1°C
- Flow rate: 8 litres per minute
- Integrated with NeuroEquilibrium VNG

NeuroEquilibrium™ provides systems for comprehensive management of vertigo and balance disorders by breakthrough innovations. Our strategy is to set up specialized advanced vertigo and balance clinics across the world in collaboration with leading ENT and neurology clinics and multi-speciality hospitals through:

- Neurotological evaluation with our suite of world class diagnostic equipment
- A team of neurotological specialists.

We are India's first robust diagnostic platform helping doctors diagnose and treat vertigo and balance disorder across the world.

Let's **get together** to
restore **balance**



**Neuro
Equilibrium**
Diagnostic Systems Pvt Ltd.

140, Sri Gopal Nagar,
Gopalpura Bypass Road,
Jaipur - 302015
India

info@neuroequilibrium.in
www.neuroequilibrium.in



An ISO 9001 and ISO 13485 certified company