


The logo consists of three overlapping loops: a blue one at the top, a black one on the left, and a green one on the right. A registered trademark symbol (®) is located to the upper right of the loops.

Neuro Equilibrium

**Cutting-edge technology for
diagnosis and treatment of
vertigo & balance disorders**

SUBJECTIVE VISUAL VERTICAL



For detection of abnormal subjective tilt

Subjective Visual Vertical (SVV) is an investigation to evaluate the otolith system, which is responsible for perception of verticality. Static and dynamic SVV are important to:

- Assess otolithic disorders
- Assess chronic dizziness
- Differentiate peripheral from central vestibular disorders
- Decide side of peripheral vestibular insult during the acute stage
- Diagnose compensated vestibular disorders
- Assess effect of rehabilitation

Product Specifications

- Tubular vision goggle with 9 axis IMU
- Real time data analysis
- Fully integrated and computerized
- Simple to use, no moving parts

Tests Conducted

- Static SVV, head vertical
- Static SVV, head inclined at 30 degrees
- Dynamic moving background, head vertical

NeuroEquilibrium™ provides systems for comprehensive management of vertigo and balance disorders by breakthrough innovations. Our strategy is to set up specialized advanced vertigo and balance clinics across the world in collaboration with leading ENT and neurology clinics and multi-speciality hospitals through:

- Neurotological evaluation with our suite of world class diagnostic equipment
- A team of neurotological specialists.

We are India's first robust diagnostic platform helping doctors diagnose and treat vertigo and balance disorder across the world.

Let's **get together** to
restore **balance**



**Neuro
Equilibrium**
Diagnostic Systems Pvt Ltd.

140, Sri Gopal Nagar,
Gopalpura Bypass Road,
Jaipur - 302015
India

info@neuroequilibrium.in
www.neuroequilibrium.in



An ISO 9001 and ISO 13485 certified company